

Test Kit, Boiler Water Model AL94 MG-L (Drop Count Titration)

Components for Kit 222403

222403	222403-EU*	222403-TH*	Description
Part Number	Part Number	Part Number	
43800	43800	43800	TUBE, MEASURING, PLASTIC 5.83ML
173000	173000	173000	TUBE, VIEWING 5-10ML-LP W/LOGO PK/6
173106	173106	173106	STOPPER, BLUE RUBEBR PK/6
173200	173200	173200	COLOR COMPARATOR BOX
232706	232706	232706	BOTTLE, MIXING MARKED PK/6
4660500	4660500	4660500	CASE ASSY, HC-5**
990100	990100	990100	COLOR DISC 2-PIECE, WIDE RANGE pH
9262400	9262400	9262400	INDOOR LIGHT, COLOR DISC 2PC, IRON-FERROVER MD-RNG
9263800	9263800	9263800	COLOR DISC 2PC, IRN-FERROVER MR OL
94299	94299-EU	94299-TH	PHENOLPHTHALEIN PP PK/100
94399	94399	94399-TH	BROMCRESOL GREEN -METHYL RED PP PK/100
2349732	2349732	2349732-TH	SULFURIC ACID, 0.035N 100ML MDB
104399	104399	104399-TH	CHLORIDE 2 IND PWD PLWS PK/100
42432	42432	42432-TH	BUFFER SOLN, HDNS 1 100ML MDB
42532	42532-EU	42532-TH	HARDNESS 2 TEST SOLN 100ML MDB
92799	92799	92799-TH	FERROVER PWD PLWS 5ML PK/100
2329332	2329332	2329332-TH	WIDE RANGE-4 PH IND, 100ML MDB
2349832	2349832	2349832-TH	SILVER NITRATE TITRANT 0.0493N 100ML
2349932	2349932	2349932-TH	EDTA STD SOLN, 0.035N 100ML MDB
DOC326.97.00087	DOC326.97.00087	DOC326.97.00087	Downloadable from Hach.com

*EU=European Union, TH=Thailand

**Case only.